


## PRINCIPLE:



A laser source is spatially filtered via a microscope objective and a pin hole. The pin hole is located at the focal point of a collimating lens. The collimated beam encounters a slightly wedged glass plate. The surface adjacent to the collimating lens is of good optical quality. However, the next surface is of exceptional optical quality,  $\lambda/10$  peak to valley (PV) or better. This is the reference surface. A part of the collimated beam is reflected by this surface. And a part of the collimated beam continues on to interrogate the component being tested.

The return beam contains information on aberration introduced by the test surface. The two wave fronts recombine inside the interferometer.

The cube beam-splitter diverts the combined beam towards a recording medium, either film or a CCD. An intermediate lens forms an image of the test surface onto the recording plane.

The images so formed are called interference fringes. The shape of the light bands (fringes) represents the true profile of the surface under test.

With the help of Fringe analysis software the flatness of the test part could be quantified and report may be generated.

### Applicaton

Ceramic Seals	Mechanical seals	Pump valves
Clutch plates	Gear Wheels (Flat-honed)	Auto Parts
Optical Flats	Optical Prisms	Semiconductor wafers

### About us

Custom built manufacturer of Plano Optics: Prisms, Windows, Flats & Mirrors.

Pioneer in Prism production for high tolerances by Optical Contacting Technique in India

Capable to produce 5 mm to 300 mm Optics within 2 Arc seconds tolerances & upto  $\lambda/10$  flatness.

In-house capabilities to offer coated Optics – Antireflection Coatings & Aluminum Reflective Coatings.

1991- Established Optics manufacturing facility.

1994- Developed Laser Fizeau Interferometer

1996- Established Optical Contacting technique for production of Prisms

2002- Developed Autocollimator/Angle Dekkor for Machine Tool Industry

2004- Manufacture of Monochromatic Light Source

2005- Manufacture of Michelson Interferometer

2005- Introduced double-side polishing machine for achieving < 1 sec parallelism.

2005- Developed Interference Autocollimator

2006- Participated in Hannover Fair, Germany.

2006- Incorporated Software evaluation to certify Flatness of Optical Flats

2007- Developed Alignment Autocollimator

2007- Participated in MTA Exhibition, Singapore

2008- Designed & Developed Optical Head for Alignment application



## PRISMS INDIA PRIVATE LIMITED

C-92, Pipdic Industrial Estate, Mettupalayam, Pondicherry – 605 009

Tel: 0413-2277059, 2271289, 4204060, Fax: 4204060, 2250328

E-mail: prisms@vsnl.com, prisms@airtelmail.in Website: www.prismsindia.com

## Flatness Measurement

# FIZEAU INTERFEROMETER

Technical mass products with flat surfaces specified to better than  $0.1\mu\text{m}$  over a few cm are produced by millions per month.

The production processes are lapping & polishing – and, increasingly more often, micro-machining; 90% have diameters under 50mm.

ISO 9000 requires tight inspection often 100%

The most common tool for this is The Interferometer

An esteemed value to engineers as does the stethoscope for the medical practitioner.

Fizeau Interferometer  
Model FIM-100



Fizeau Autocollimator  
Model FAC-30, 60




### Specifications

Orientation: Vertical  
Clear Aperture: 95 mm  
He-Ne Laser (633 nm): 1.0 mW  
Reference Flat:  $\lambda/20$   
Work Table (Tilt Table) Dimension: 6" x 6"  
1/3" CCD Camera  
Spot alignment & Direct Fringe viewing  
Dimensions: 350 (L) x 400 (W) x 670 (H)  
Weight: 30 Kgs

Orientation: Vertical  
Clear Aperture: 55 mm  
Diode Laser (650 nm): < 3mW  
Reference Flat:  $\lambda/10$   
Work Table (Tilt Table) Dimension: 4" x 4"  
1/3" CCD Camera  
Direct Fringe viewing  
Dimensions: 300 (L) x 250 (W) x 550 (H)  
Weight: 20 Kgs

### Optional Accessories

Static Quick Fringe Analysis Software & Frame Grabber



## PRISMS INDIA PRIVATE LIMITED